

Investigation of metabolic changes due to recombinant protein production in *E. coli* using 2D fluorescence spectroscopy and principal component analysis

Garima Jain, Mickie Takagie, Ursula Rinas, Guhan Jayaraman, Bernd Hitzmann

Institut für Technische Chemie, Callinstr. 3, 30167 Hanover, Germany
Dept. of Biotechnology, Indian Institute of Technology, 600036 Madras, India
Helmholtz Centre of Infection Research, 38124 Braunschweig, Germany

INTRODUCTION: 2-Dimensional fluorescence spectroscopy can be used as an on-line tool for continuous monitoring of bioprocesses. It monitors the metabolic state of cells on-line and gives useful information about various key molecules, e.g. FAD/FMN, NADH, pyridoxine, tryptophan, tyrosine and phenylalanine. This study gives a comparison in the metabolic state of *Escherichia coli* [BL21(DE3)+pET-29c-hbFGF] with and without recombinant protein (human basic fibroblast growth factor) production in batch cultivations using principal component analysis (PCA).




Figure 1: Showing an illustration of principal component analysis

Principal component analysis of induced and uninduced cultures:




Figure 2

Figure 3

Figs. 2 and 3: Showing dry cell weight (g/L) (DCW), specific carbon dioxide production rate (mM/g/h) (Sp CPR), principal component 1 (PC 1) and principal component 2 (PC 2) values for uninduced and induced batch cultivations, respectively.

OBSERVATIONS AND CONCLUSIONS :

- PC 1 component correlates very well with biomass formation.
- PC 2 correlates with the specific carbon dioxide production rate in the batch cultivation with recombinant protein production. This shows that PC 2 contains information about major metabolic changes after induction.




Figure 4 Induced batch

Figure 5 Uninduced batch

Fig. 4 and 5: Showing surface plots of fluorescence spectrum (at points with same biomass) at the time of induction (DCW = 0.8 g/L) and spectrum at the end of cultivations (DCW = 3.0 g/L) in induced and uninduced batch cultivations, respectively.

OBSERVATIONS AND CONCLUSIONS :

- At the end of the cultivation, the amount of FMN/FAD and NADH is relatively higher in batch cultivation with induction as compared to batch cultivation without induction.
- Amount of FMN/FAD and NADH are affected due to changes in the respiratory activity upon induction.
- This is indicated by the correlation of PC 2 with the specific carbon dioxide production rate as shown in Fig. 3.

SUMMARY: 2-D fluorescence spectroscopy can give useful insights in metabolic changes in *E. coli* in response to recombinant protein production. PC 2 data alignment with sp. CPR (only in induced batches) indicates that PCA contains information related to changes in the respiratory activity. These changes are mainly observed in the FMN/FAD and NADH region of the spectrum.

FUTURE PLANS:

- To make partial least square-regression models for prediction of various off-line measurements e.g. biomass, glucose, acetate.
- To extend this analysis for fed-batch cultivations.

ACKNOWLEDGEMENT: G. Jain would like to thank the DAAD for a fellowship.